


CATALOG NO.	NOMINAL CAPACITY, (mL)							MATERIALS	
								VIAL	CLOSURE
006250	5	12	18	28	50	75	160	LDPE	LDPE

SIZE CODE.....:	-0005	-0012	-0018	-9028				
SIZE (NOM VOL mL) ..:	5	12	18	28				
NOM. BRIM CAP. (mL):	6.2	12.7	19.7	29.8				
NOM. WEIGHT W/CLOSURE (gms).....:	2.3	3	4	7				
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
	0.74±.02	18.8±0.5	0.88±.02	22.4±0.5	0.88±.02	22.4±0.5	0.88±.02	22.4±0.5
	0.97±.03	24.6±0.8	1.39±.05	35.3±1.3	2.09±.05	53.1±1.3	3.35±.06	85.1±1.5
	0.93±.03	23.6±0.8	1.36±.05	34.5±1.3	2.05±.05	52.1±1.3	3.31±.06	84.1±1.5
	0.93±.03	23.6±0.8	1.07±.04	27.2±1.0	1.07±.04	27.2±1.0	1.07±.04	27.2±1.0
	0.79±.02	20.1±0.5	0.92±.03	23.4±0.8	0.91±.03	23.1±0.8	0.89±.03	22.6±0.8
	0.035	0.89	0.035	0.89	0.035	0.89	0.040	1.02

SIZE CODE.....:	-9050	-9075	-0160			
SIZE (NOM VOL mL) ..:	50	75	160			
NOM. BRIM CAP. (mL):	52	80	168			
NOM. WEIGHT W/CLOSURE (gms).....:	10	14	23			
	INCHES	MM	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	1.26±.03	32.0±0.8	1.26±.03	32.0±0.8	1.88±.04	47.8±1.0
B (HT. W/CLOSURE)...	2.79±.05	70.9±1.3	4.24±.06	107.7±1.5	3.92±.06	99.6±1.5
C (HT. W/O CLOSURE):	2.74±.05	69.6±1.3	4.19±.06	106.4±1.5	3.87±.06	98.3±1.5
D (CLOSURE O.D.)...	1.48±.04	37.6±1.0	1.48±.04	37.6±1.0	2.12±.04	53.8±1.0
E (BASE O.D.).....:	1.29±.04	32.8±1.0	1.27±.04	32.3±1.0	1.91±.04	48.5±1.0
W.T. (NOM WALL THK):	0.045	1.14	0.045	1.14	0.050	1.27

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Ponorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-6250-XXXX-1	REV.	0795
CHK'D	PVC		7/13/95			
DESCRIPTION: SAMPLE VIALS WITH CLOSURES						