


CATALOG NO.....:		6420-0500		6420-0750	
MATERIALS	SPIGOT....:	LDPE		LDPE	
	O-RING....:	TEFLON TFE		TEFLON TFE	
THREAD SIZE		1/2 NPT		3/4 NPT	
NOM. WT. (gms).....:		66		72	
		INCH	MM	INCH	MM
A (OVERALL LENGTH)..:		4.2	107	4.6	117
B (OVERALL HEIGHT)..:		3.1	79	3.1	79
C (HANDLE DIA).....:		2.3	58	2.3	58



SPECS ON BACK

NOT TO SCALE

- NOTES:
- 1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 - 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC	APP'D PVC	PSS 00-6420-0XX0	REV. 0601
CHK'D PVC	6/8/01		
DESCRIPTION: NEEDLE-TYPE TANK SPIGOTS			