


| CATALOG NO. | MATERIALS |          |
|-------------|-----------|----------|
| 6700        | SHIELD    | FEET     |
|             | ACRYLIC   | NEOPRENE |

| SIZE CODE.....         | -1812  |     | -2418  |     |
|------------------------|--------|-----|--------|-----|
|                        | INCHES | MM  | INCHES | MM  |
| A (WIDTH).....         | 12     | 305 | 18     | 457 |
| B (DEPTH).....         | 12     | 305 | 12     | 305 |
| C (HEIGHT W/FEET)..... | 18.1   | 460 | 24.1   | 612 |
| D.....                 | 8      | 203 | 10     | 254 |
| E (THICKNESS).....     | 0.38   | 9.6 | 0.38   | 9.6 |

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.
1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



LOGO, TRADENAME AND CATALOG NUMBER PRINT

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|       |     |       |        |                  |              |
|-------|-----|-------|--------|------------------|--------------|
| MADE  | PB  | APP'D | PVC    | PSS 00-6700-XXXX | REV.<br>0402 |
| CHK'D | PVC |       | 4/2/02 |                  |              |

DESCRIPTION:  
BENCHTOP BETA SHIELDS