


CATALOG NO.	NOMINAL SIZE (mm)	MATERIAL
007130	376 x 241 x 114	WHITE HDPE
* 007134	376 x 241 x 114	WHITE PP

* CAT. NO. 007134 IS AUTOCLAVABLE

SIZE CODE.....:	-0010	
NOM. WEIGHT (gms)....:	440	
	INCH	MM
A (OVERALL).....:	14.8	376
B (OVERALL).....:	9.5	241
C (OVERALL).....:	4.5	114
D (TOP OF HANDLE)....:	6.1	155
X (COMPARTMENT BASE):	13.3	337
Y (COMPARTMENT BASE):	3.3	83
Z (COMPARTMENT HT.)..:	4.4	111



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE IN CHINA	JGJ	APP'D PVC	6/19/95	PSS 00-713X-XXXX-1	REV. 0695
------------------	-----	--------------	---------	--------------------	--------------

DESCRIPTION:
UTILITY CARRIERS