


CATALOG NUMBER.....:	14100-0002		14100-0005		14100-0010		14100-0015	
NOM. CAP (GAL).....:	2		5		10		15	
MATERIAL.....:	LLDPE		LLDPE		LLDPE		LLDPE	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (LENGTH).....:	10 1/2	267	16 11/16	424	14 9/16	370	20 11/16	525
B (WIDTH).....:	10 1/2	267	12 11/16	322	14 9/16	370	14 1/2	368
C (HT. W/O COVER)....:	9 1/4	235	11 7/16	291	14 5/16	364	13 9/16	235
D (LENGTH).....:	8 3/16	208	14 1/2	368	12 1/4	311	18 7/16	344
E (WIDTH).....:	8 3/16	208	10 3/8	264	12 1/4	311	12 1/2	318
F (LENGTH COVER).....:	11 3/16	284	17 5/16	440	15 1/16	383	21 5/16	541
G (WIDTH COVER).....:	11 3/16	284	13 5/16	338	15 1/16	383	15 1/4	387
H (HT. COVER).....:	1 1/4	32	1 1/4	32	1 1/4	32	1 1/4	32
W.T. (NOM. WALL THK.)	5/32	4	5/32	4	5/32	4	5/32	4

CATALOG NUMBER.....:	14100-0020		14100-0040		14100-0045		14100-0065	
NOM. CAP (GAL).....:	20		40		45		65	
MATERIAL.....:	LLDPE		LLDPE		LLDPE		LLDPE	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (LENGTH).....:	20 11/16	525	26 5/8	676	26 5/8	676	26 11/16	678
B (WIDTH).....:	14 5/8	371	14 1/2	368	20 11/16	525	20 9/16	522
C (HT. W/O COVER)....:	19 13/16	503	13 1/4	337	19 9/16	483	25 3/16	640
D (LENGTH).....:	18 7/16	468	24 3/8	619	24 3/16	614	24 7/16	621
E (WIDTH).....:	12 1/2	318	12 1/2	318	18 7/16	468	18 5/16	465
F (LENGTH COVER).....:	21 5/16	541	27 1/8	689	27 7/16	697	27 7/16	697
G (WIDTH COVER).....:	15 1/4	387	15 1/8	384	21 1/8	537	21 1/8	537
H (HT. COVER).....:	1 1/4	32	1 1/4	32	1 1/4	32	1 1/4	32
W.T. (NOM. WALL THK.)	5/32	4	5/32	4	5/32	4	5/32	4

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	PWC	PSS 14100-00XX	REV. 0314
CK'D	KOW	3/27/14			

DESCRIPTION: TANKS WITH COVERS