



CATALOG NUMBER.....:	312160-0200	312160-0240	312160-0384	312160-0530				
MATERIAL.....:	PVC, CLEAR	PVC, CLEAR	PVC, CLEAR	PVC, CLEAR				
T (THICKNESS)(mil)::	3*		3*		3*		4*	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
W (FLAT WIDTH).....:	1.57 ±.04	40 ±1	1.81 ±.04	46 ±1	2.72 ±.04	69 ±1	4.10 ±.04	104 ±1
H (HEIGHT).....:	1.1 ±.04	27 ±1	1.2 ±.04	30 ±1	1.9 ±.04	47 ±1	1.6 ±.04	40 ±1

\* THICKNESS TOLERANCE = ±10%



SLEEVE ILLUSTRATION  
NO SCALE



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PSS 312160-0XXX	REV.
CK'D	CG	BR 5/31/06		0506

DESCRIPTION:  
PVC SHRINKBAND