

## Thermo Scientific Abgene® Blank and Alphanumeric Storage Tubes

Thermo Scientific Abgene Blank and Alphanumeric Storage Tubes enable secure sample storage, while allowing individual access, sealing and retrieval of samples. Storage tubes can be picked and placed for reformatting of sample selections and eliminate the need to thaw an entire storage plate for only a few desired samples.



**Highest Quality.** All Thermo Scientific Abgene blank and alphanumeric storage tubes are manufactured from 100% virgin, medical-grade polypropylene. Storage tubes are leak tested to ensure the integrity and security of samples and with a V-bottom design, our alphanumeric storage tubes facilitate sample recovery.

**Enhanced Rack Design.** Unlike traditional storage tubes or blocks, our blank and alphanumeric storage tubes are available in our specially-designed, stackable, 96-format, microplate-footprint racks to save precious space in your freezer, refrigerator, compound store or on the bench top. In addition, the rack design includes grid referenced lids for simple identification of tube location.

**Superior Sample Protection.** Several sealing options are available in sterile and non-sterile, solid or pre-split cap designs, both pre-installed or in bulk for fitting at the bench. Pierceable cap designs, such as

Thermo Scientific Abgene SmartCap Mats and the Thermo Scientific Matrix® SepraSeal® capping system, offer the flexibility to seal and subsequently access sealed individual storage tubes. For applications involving repetitive sample access, maintain sample integrity with Thermo Scientific Abgene Multisip™ Plugs and Thermo Scientific Matrix DuraSeal® caps. With the Thermo Scientific Multisip Applicator or Thermo Scientific SuperSealer ensure sealing is secure and consistent every time. Additionally, Thermo Scientific Abgene Twist-Lock™ caps are available for most volumes.

**Knowledgeable Support.** Our team has a strong reputation for supporting both equipment and consumables, and can demonstrate innovative ways to store samples in your laboratory using our blank and alphanumeric storage tubes.

## Thermo Scientific Abgene Blank and Alphanumeric Storage Tubes

Tube Material	Polypropylene
Rack Material	Polypropylene
Autoclavable	Yes

### Ordering Guide

Item No.	Description	Qty
AB-1417	Abgene 0.3 ml Blank Tubes	10 racks of 96 tubes/case
AB-1418	Abgene 0.3 ml Blank Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1264	Abgene 0.65 ml Blank Tubes, Bulk	1000 tubes/case
AB-1256	Abgene 0.65 ml Blank Tubes	10 racks of 96 tubes/case
AB-1257	Abgene 0.65 ml Blank Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1429	Empty Racks for Abgene 0.65 ml Blank Tubes, Purple	10 racks/case
AB-0672	Abgene 1.2 ml Blank Tubes, Bulk	1000 tubes/case
AB-0595	Abgene 1.2 ml Blank Tubes	10 racks of 96 tubes/case
AB-0970	Abgene 1.2 ml Blank Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1432	Empty Racks for Abgene 1.2 ml Blank Tubes, Purple	10 racks/case
AB-1407	Abgene 0.3 ml Alphanumeric Tubes	10 racks of 96 tubes/case
AB-1408	Abgene 0.3 ml Alphanumeric Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1259	Abgene 0.65 ml Alphanumeric Tubes	10 racks of 96 tubes/case
AB-1260	Abgene 0.65 ml Alphanumeric Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1430	Empty Racks for Abgene 0.65 ml Alphanumeric Tubes, Red	10 racks/case
AB-0977	Abgene 1.2 ml Alphanumeric Tubes	10 racks of 96 tubes/case
AB-0979	Abgene 1.2 ml Alphanumeric Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1228	Abgene 1.2 ml 2D Barcoded and Alphanumeric Tubes	10 racks of 96 tubes/case
AB-1229	Abgene 1.2 ml 2D Barcoded and Alphanumeric Tubes w/ Multisip Plugs installed	10 racks of 96 tubes/case
AB-1433	Empty Racks for Abgene 1.2 ml Alphanumeric Tubes, Red	10 racks/case
AB-1461	Abgene 0.5 ml Blank Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1459	Abgene 0.5 ml Blank Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1263	Abgene 1 ml Blank Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1262	Abgene 1 ml Blank Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1316	Abgene 0.5 ml Alphanumeric Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1315	Abgene 0.5 ml Alphanumeric Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1261	Abgene 1 ml Alphanumeric Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1258	Abgene 1 ml Alphanumeric Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1250	Abgene Twist-Lock Caps, Bulk	960 caps/case
AB-1317	Abgene Colored Twist-Lock Caps, Bulk, 240 each of 4 colors	960 caps/case
AB-1255	Abgene Twist-Lock Tool	Each
AB-1312	Abgene Twist-Lock Cap-It Removal Tool	Each

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: [matrix.info@thermofisher.com](mailto:matrix.info@thermofisher.com)  
 Europe: Tel: +44 (0) 161 486 2110 | email: [matrix.eu.info@thermofisher.com](mailto:matrix.eu.info@thermofisher.com)  
 Asia: email: [matrix.ap.info@thermofisher.com](mailto:matrix.ap.info@thermofisher.com)

[www.thermo.com/matrix](http://www.thermo.com/matrix)

**Thermo**  
SCIENTIFIC