












Thermo Scientific™ Samco™ Exact Volume Transfer Pipette

- ▶ Polyethylene, non-sterile
- ▶ Useful for small samples or dispensing reagents
- ▶ Ideal for home & POC (Point of Care) testing
- ▶ Useful in clinical laboratories
- ▶ OEM applications
- ▶ Custom design options available to meet your specifications

Cat. No. **Dispense volume (µl)** **Stem Diameter (mm)** **Length (cm)** *Exact volume pipette lengths indicated below are close to actual size*

783NL	20	2.0	5.28	
784NL	40	2.3	5.92	
930NL	50	2.5	5.92	
785NL	60	2.5	6.48	
786NL	80	2.5	7.49	
787NL	100	2.5	8.30	
788NL	150	3.0	7.95	
941NL	200	3.0	8.89	
191NL	250	3.0	10.19	
966NL	300	3.6	9.50	
967NL	400	3.8	10.59	

500 pipets per pack in all sizes. 5000 pipets per case in all sizes.

thermoscientific.com

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Scientific Customer Service:
Tel: 1-800-995-2787 / Fax: 1-858-546-1642
E-mail: Info.sandiego@thermofisher.com

BRCCSSAMCOEPIIP 0615

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand