

Specialized Chemical Services

Your requirement is our goal—
turn to a trusted chemicals partner



Looking for a chemistry solution?



When products you need are not readily available off the shelf from a catalog, you can save time and resources by turning to the company you can trust for your custom chemical requirements.

Can I get this product in a different pack size from what's in your catalog?

Yes, we can adapt our packaging to your exact needs

Can you prepare the solvent blend I need?

Yes, we can help you save time and effort by providing you with a custom-made mixture

Can I include another specification that's important to me?

Yes, we can adapt our specifications by adding the parameter required

I need 397 lbs of a product; is this possible?

Yes, you can drive efficiency with our dedicated packaging size

Can I add any additional details to the label?

Yes, we can customize our label to meet your needs and enable batch-to-batch consistency

Optimize your bench-to-production workflow

How can we help?

We customize our products to meet your requirements, offering you an ideal solution. We distill, blend, test, package, and deliver chemicals tailored to your needs.

Discover the extensive benefits of Fisher Chemical™ customized products:

Semibulk and bulk chemical services

Benefits:

- Simplify scale-up cost-effectively
- Drive efficiency

Custom-made chemicals

Benefits:

- Help save time and effort, while reducing waste
- Enhance reliability
- Enable consistent quality
- Provide convenient handling, storage, and preservation of product integrity

Returnable drum delivery system

Benefits:

- Helps reduce cost and waste
- Helps reduce exposure to hazardous materials
- Provides a recycling solution



Semibulk and bulk chemical services


Looking to scale up?

Looking to simplify the process by using high-quality chemicals?

Looking to adapt packaging to scale?

With our extensive portfolio of analytical chemicals, we offer one of the most comprehensive, cost-effective solutions available in the market for all your production chemical needs. Whether you are engaged in initial research or scaling up to full production, we can customize products exclusively for you.

Drive efficiency by moving to semibulk and bulk solutions for your solvents and inorganics.



Custom-made chemicals



Standard catalog chemicals customized to your requirements

Tailored solvents and custom blends

Let us prepare chemical blends tailored to your needs, to help you save time, reduce waste, and achieve consistency.

Customized specification and testing

Optimize and enhance the reliability of your process by taking advantage of our extensive analytical services to meet your specifications, which will also be included in your certificate of analysis.

Custom packaging and labeling

Help improve handling, storage, and product integrity by receiving the exact amount of material you want in the exact packaging you want.

FisherPak Solvent Delivery System

The Fisher Chemical™ FisherPak™ Solvent Delivery System is an ideal choice for handling large quantities of high-purity solvents. Almost 200 lbs of glass and cardboard are eliminated by each 200 L FisherPak system (available in 19 L, 50 L, 200 L, and 1,350 L sizes). From lab to production scale, the FisherPak solution delivers solvents with enhanced safety and improved productivity.

The FisherPak Solvent Delivery System protects solvent purity throughout shipment, delivery, installation, and dispensing, offering:

- Pressurized Type 316 stainless steel electropolished containers with superior corrosion-resistance properties
- Stainless steel cover and tamper-evident closure
- An enclosed system that includes chemical-dedicated pumps and lines to minimize contamination of solvent during dispensing
- All containers dedicated to a specific customer and solvent at all times

Improves productivity

The FisherPak Solvent Delivery System helps to increase lab efficiency while reducing lab operating costs by:

- Reducing solvent testing
- Eliminating bottle rinsing
- Cutting disposal costs
- Providing a wide range of returnable containers to meet your needs
- Offering extra testing parameters, blending, and customized capabilities for mechanical, packaging, and logistics solutions

Enhances safety

The FisherPak system enables safe practices through:

- A unique combination of mechanical and manual controls to help protect the lab and the environment
- Its bottle-free, closed system that eliminates glass breakage and helps reduce the risk of spills and exposure to vapors

- Unique horizontal vapor venting for superior fire safety
- Secondary manual shut-off valves, in addition to isolation valves, incorporated into the quick-connect hardware
- Color-coded mechanical lockout connectors that help prevent improper hookups
- Compliance with Department of Transportation and NFPA 30 standards
- ASME (American Society of Mechanical Engineers) certification of returnable containers

Promotes environmentally friendly practices


The FisherPak system can help to reduce the amount of solid waste generated in your laboratory, thereby reducing environmental impact.


- A reusable 200 L FisherPak container replaces fifty 4-liter bottles and packing material
- Flammable or toxic solvent liquids and vapors are not released into the surrounding air
- Empty FisherPak containers are cleaned and refilled by Thermo Fisher Scientific—there are no more bottles or associated risks



How does the cycle work?

- 1 Project discovery—our team of specialists will assist in selecting the right drum type and accessories, and in calculating the necessary fleet size to meet your specific solvent and delivery requirements
- 2 Process initiation—rotation process defined and quote generated
- 3 First order placed
- 4 Shipment of material
- 5 Use material until container is empty; contact us to arrange container collection
- 6 Container cleaned and validated
- 7 Container refilled and released for shipment
- 8 Project optimization—usage and supply chain are monitored and adjusted to reflect changes in demand and/or logistics





The chemicals you need—tested,
packaged, and delivered to meet your
exact requirements.

Let our experienced team make your
life easier.

Request a quote at
fishersci.com/bulkcustomchemicals



In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: fishersci.com

In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: fishersci.ca



For Research Use or Further Manufacturing. Not for diagnostic use or direct administration into humans or animals.

© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **BN040317CH83 1118**