

Rotors and Carriers Ordering Information

SpeedVac Models		DNA 120	ISS110	SPD111V	SC210A	Rotor Part Number		Rotor System	
		DNA 1200P	SPD1010	SPD121P	SPD2010				
Microtubes	No. of Tubes	Small Capacity	Medium Capacity		Large Capacity				
0.5ml	48	X					RD48		Traditional
0.5ml	72	X					RD72		Traditional
1.5ml	24	X					RD24		Traditional
1.5ml	36	X					RD36		Traditional
1.5ml	40		X	X			RH40-11		Traditional
1.5ml	50				PR*	PR*	RBA25-11-39		Set of 2 blocks for PRO
1.5ml	64		X	X			RH64-11		Traditional Rotor System
1.5ml	100				PR*	PR*	RBA25-11-39 (qty 2)		Set of 2 blocks for PRO
1.5ml	120		X	X			RH120-11		Traditional
1.5ml	200				X	with RAD	RH200-12		Traditional
Glass and Plastic Tubes									
6 x 32 mm	252				PR*	PR*	RBA126-6-32		Set of 2 blocks for PRO
6 x 32 mm	504				PR*	PR*	RBA126-6-32 (qty 2)		Set of 2 blocks for PRO
6 x 50 mm	100		X	X			RH100-6		Traditional
8 x 29 mm	100		X	X			RH100-8		Traditional
12 x 75 mm	40		X	X			RH40-12		Traditional
12 x 75 mm	72		X	X			RH72-12		Traditional
12 x 75 mm	200				X	X (requires RAD200)	RH200-12		Traditional
13 x 100 mm	32		X	X			RH32-13		Traditional
13 x 100 mm	70				PR*	PR*	RBA35-13-100		Set of 2 blocks for PRO
13 x 100 mm	118				X	X (requires RAD200)	RH200-12		Traditional
13 x 100 mm	140				PR*	PR*	RBA35-13-100 (qty 2)		Set of 2 blocks for PRO
16 x 100 mm	48				PR*	PR*	RBA24-16-100		Set of 2 blocks for PRO
16 x 100 mm	98				PR*	PR*	RBA24-16-100 (qty 2)		Set of 2 blocks for PRO
16 x 125 mm	48				X	X (requires RAD200)	RH48-18-125		Traditional
16 x 125 mm	30				PR*	PR*	RBA15-16-125		Set of 2 blocks for PRO
17 x 60 mm	48				PR*	PR*	RBA24-17-60		Set of 2 blocks for PRO
17 x 60 mm	96				PR*	PR*	RBA24-17-60 (qty 2)		Set of 2 blocks for PRO
16 x 125 mm	60				PR*	PR*	RBA15-16-125 (qty 2)		Set of 2 blocks for PRO
17.5 x 102 mm, 15 ml Corex®	8		X	X			RH8-17.5		Traditional
18x100,17x95, 16x100mm	8		X	X			RH8-18		Traditional
18 x 125 mm	48				X	X (requires RAD200)	RH48-18-125		Traditional
18 x 150 mm, 16 x 125 mm	32				X	X (requires RAD200)	RH32-18-150		Traditional
28 x 150 mm	12				X		RH12-29		Traditional

Rotors and Carriers Ordering Information Continued

SpeedVac Models		DNA 120		SPD111V		SC210A		
		DNA 1200P	ISS110	SPD121P	SPD1010	SPD131DDA	SPD2010	SC250EXP
Centrifuge Tubes	No. of Tubes	Small Capacity	Medium Capacity		Large Capacity		Rotor Part Number	Rotor System
15ml conical	10		X	X	X		RH10-15	PR*
15ml conical	30				X	X (requires RAD200)	RH60-17-100	Traditional
50 ml conical	6		X	X			RH6-50	Traditional
50 ml conical	12				X		RH12-29	Traditional
Flasks								
100 ml pear-shaped flask	4		X	X			RH4-100	Traditional
100 ml pear-shaped flask	8				X		RH8-200	Traditional
Vials								
12 x 32 mm	60		X	X			RH60-12-40	Traditional
1 dram (15 x 45 mm), 4 ml	24		X	X			RH24-15	Traditional
1 dram (15 x 45 mm), 4 ml	192				X	X (requires RAD200)	RH192-15	Traditional
12 x 40 mm	60		X	X			RH60-12-40	Traditional
20 x 47 mm	12		X	X			RH12-20	Traditional
20 x 60 mm	12		X	X			RH12-20	Traditional
18 x 52 mm mini-scintillation	24		X	X			RH24-18	Traditional
28 x 58 mm, 20 ml scintillation	24				PR*	PR*	RBA12-28-58	Set of 2 blocks for PRO
28 x 58 mm, 20 ml scintillation	48				PR*	PR*	RBA12-28-58 (qty 2)**	Set of 2 blocks for PRO
28 x 60 mm scintillation	12		X	X			RH12-28	Traditional
28 x 60 mm scintillation	50				X	X (requires RAD200)	RH50-28-60	Traditional
19 x 65 mm flat bottom	48				PR*	PR*	RBA24-19-65	Set of 2 blocks for PRO
19 x 65 mm flat bottom	98				PR*	PR*	RBA24-19-65 (qty 2)**	Set of 2 blocks for PRO
Microwell Plates								
Shallow well plates	2	X					RD2MP	Traditional
Shallow well plates	12				X		MPTR12-210	Traditional
Shallow well plates	4				PR*	PR*	UPC-1	Set of 4 carriers for PRO
Shallow well plates	6			X			RHSW6MP	Traditional
Shallow well plates	2		X	X			RH2MP	Traditional
Shallow well plates	8				PR*	PR*	UPC-2	Set of 4 carriers for PRO
Shallow well plates	20				PR*	PR*	UPC-5	Set of 4 carriers for PRO
Deep well plates	2			X			RHDW2MP	Traditional
Deep well plates	4				PR*	PR*	UPC-1	Set of 4 carriers for PRO
Deep well plates	8				PR*	PR*	UPC-2	Set of 4 carriers for PRO
Deep well plates	8				X		MPTR8-210	Traditional

* For use with the PRO rotor system

** To achieve the tube capacity listed, please order quantity 2 of the rotor blocks