

Thermo Scientific Smart-Vue Wireless Monitoring Solution

Technical Data Sheet: Telephone Modem



Application

USB Telephone modem

General

- US Robotics 56K* USB Modem
- Phone cord with RJ11 terminal connector
- USB extension cable
- Velcro* Strap for cable management
- CD-ROM with drivers (P/N:SVU601-LSB)
- Package Dimensions/Weight: 260mmX160mmX 0mm, (10.24"X6.30"X2.36"). Weight 0.3
- Module Holder Mounting Options: Velcro*, magnetic, or fixed
- Product dimensions: 3.88 x 1.5 x .88 in (9.86 x 3.81 x 2.24cm)
- Weight: 0.26 lb (0.12 kg)
- Warranty: 1 year on all hardware (see warranty terms & conditions for details on coverage)

Standards and Regulations

- Regulatory: FCC approved, IC approved, UL/CUL listed, TBR-21

Environmental Conditions

- Operating Temperature conditions: 0 to 50 degrees C.
- Non-Operating Temperature conditions: -20 and above 70 degrees C.
- Humidity: Operating Humidity conditions: 20% to 80% non-condensing
- Non-Operating Humidity conditions: 5% or above 95% non condensing

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Velcro is a registered trademark of Velcro Industries B.V. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

REV: 10/8/2012

Thermo
SCIENTIFIC