

Thermo Scientific ADA Compliant Value Lab Refrigerators and Freezers

An economical solution for safe sample storage in compliance with ADA height standards

Thermo Scientific Value lab refrigerators and freezers are particularly suited for routine undercounter sample storage, and at a height of 32", they will comply with countertops built to the height stipulations of the Americans with Disabilities Act (ADA).

Alert: These units are not intended for critical lab samples. Storage of high-value, life science materials or clinical biological may require alarm, monitoring, back-up and temperature performance features. The need to add these features should be determined by your industrial hygienist or safety officer. These units are not suitable for corrosive environments.

- Factory tested prior to shipment
- NRTL listed
- Environmentally-friendly, non CFC refrigerants
- Adjustable temperature control
- Easy-to-clean, durable finish
- Manual defrost

Warranty: 13 months parts and labor

Certifications: UL

**Available in North America and select Latin America countries,
consult your sales representative for availability.**



05LRAETSA open

Thermo
S C I E N T I F I C

Thermo Scientific ADA Compliant Value Lab Refrigerators (Temperature Range 1° to 12°C)

Model No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Certification	Shipping Weight (lbs.)
05LRAETSA	5.5	manual	115V / 60Hz	4 adjustable	yes	17.5 x 20 x 27.75	24.25 x 23.65 x 32	UL	102

Thermo Scientific ADA Compliant Value Lab Freezers (Temperature Range -12° to -20°C)

Model No.	Capacity cu. ft.	Defrost	Voltage	Drawers	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Certification	Shipping Weight (lbs.)
04LFAETSA	3.5	manual	115V / 60Hz	3	yes	16.5 x 17.75 x 24.25	23.5 x 23.5 x 32	UL	104



04LFAETSA



05LRAETSA

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada: +1 866 984 3766 (866-9-THERMO)

Latin America: +02 28 03 20 00

DSCSADALRF 1115

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific