





Technical Data Sheet
DeviceLink™ Connect Door Switch Independent Sensor
 Thermo Fisher Scientific, Asheville, North Carolina
 Revision A

Specifications	Catalog Number	
	DLCPTDOOR001	

Application and Operating Conditions	
Application	Magnetic switch sensor to monitor door opening and closing via a DeviceLink Connect remote monitoring system
Application Environment	Non-corrosive, non-flammable, non-explosive, good air ventilation
Sensor Specifications	
Connection Interface	Digital Input (RJ45)
Type	Magnetic door switch
Operating Range	Binary (open/closed) (normally open)
Accuracy	n/a
Display Resolution	n/a
System Resolution	n/a
Cable Length	6.5' (1.98m)



© 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.