





CATALOG NUMBER.....:	5977-0013/-0313		5977-0017/-0317	
MATERIAL.....:	FILLED PP		FILLED PP	
ARRAY	8 x 12		5 x 10	
FITS TUBE SIZES (mm).....:	10 - 13		14 - 17	
NOM. WEIGHT (gm).....:	200		180	
	INCH	MM	INCH	MM
A (LENGTH).....:	7.0	178	7.4	188
B (WIDTH).....:	4.94	126	4.08	104
C (PEG HEIGHT).....:	1.75	44.5	2.00	50.8
D (OVERALL HEIGHT).....:	2.50	64	2.75	70

CATALOG NUMBER	COLOR
5977-0013	WHITE
5977-0313	BLUE
5977-0017	WHITE
5977-0317	BLUE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ISOMETRIC VIEW



5977-0017/-0317 RACK SHOWN

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-5977-XXXX	REV. 0102
CK'D	PVC		1/4/02		

DESCRIPTION: TEST TUBE PEG RACKS