


## New Improved Design!

We heard your feedback and we've made improvements to the Thermo Scientific™ ClipTip™ 1250µl tips.

The new design includes a thinner, pointier tip shape to allow more precise aspiration from the bottom of reagent vessels. We have also made the tip a bit longer.

The raw material used in the filter has been changed for filtered tips to help improve tip performance in stepping applications.

The new CT1250 tips are compatible with corresponding Thermo Scientific™ E1-ClipTip™ pipettes and have the same performance specifications as before. The product catalog numbers stay the same.



For more information, go to [thermofisher.com/cliptip](http://thermofisher.com/cliptip)

Package description	Cat. no.
ClipTip 1250, Racked, 96 tips/rack, 8 racks/package	94410810
ClipTip 1250, Racked, Sterile, 96 tips/rack, 8 racks/package	94410813
ClipTip 1250, Reload Insert, 96 tips/insert, 8 inserts/package	94410817
ClipTip 1250 Reload Insert, Sterile, 96 tips/insert, 8 inserts/package	94410818
ClipTip 1250, Racked, Filtered, Sterile, 96 tips/rack, 8 racks/package	94420813
ClipTip 1250, Reload Insert, Filtered, Sterile, 96 tips/insert, 8 inserts/package	94420818