

Snappies - Born of necessity

Collect, store and feed using a single
innovative breast milk container

A leak-resistant, airtight solution

Available in the 11 mL Colostrum Collector, 1 oz. and 2.3 oz sizes.

Our patented airtight, leak-resistant containers give lactation consultants, care givers and mothers a helping hand to make breastfeeding easier in those early, critical months.

Used in hospitals and NICUs across North America, Snappies® enable moms to collect, store and feed breast milk from the same container.



COLLECT

- **Minimize contamination:** The flip-top is conveniently hinged to the container so it will never be lost or dropped and can be easily opened and closed with one hand.
- **Direct connection:** Compatible with Ameda®, Medela® and most other standard sized pump rings.

STORE

- **Leak proof:** Our patented triple-seal design provides an airtight and leak-resistant seal for maximum shelf life in a refrigerator or freezer.
- **Ideal material:** Virgin food grade polypropylene does not contain Bisphenol A, PVC or DEHP.

FEED

- **No transfer:** The baby can be fed from the same container using a standard size bottle ring and nipple.
- **Convenient sizes:** Snappies are offered in 1oz. (35mL) and 2.3 oz. (70mL) volumes to scale with mother's milk supply.
- **Disposable and recyclable:** Making the world better for the next generation.

Features & Benefits

Leak resistant closure

Leak proof

Snappies patented triple-seal design prevents liquids from entering or exiting the container.

Flip-top hinge

Convenient

The lid locks back out of the way while in use and closes with an audible “snap” to indicate an airtight, leak-resistant seal.

Built-in label

Easy identification

Record the mother’s name, date and collection time to eliminate any potential identification problems.

Graduations

Dual graduations

Accurate, easy-to-read graduations in both oz. and mL.

Wide base

Stable

Large molded footing keeps Snappies stable while in use or when stored on the refrigerator shelf.



An all-in-one solution

Simplify the collection, storage and feeding of expressed breast milk. Snappies eliminate the need to transfer precious breast milk between containers and thus lowers the chance of contamination and risk of infection.



2.3 oz. Snappies

- Our largest Snappies container designed for when Mother's milk supply is fuller and baby's nutritional needs have increased
- Compatible with standard size bottle and pump rings

Product	Packaging	Cat No.
2.3 oz (70 mL)	200 per case – 20 packs of 10 containers	10253
2.3 oz (70 mL)	Retail pack – 1 box of 12 containers	10255
2.3 oz (70 mL)	Retail case (144 count) – 12 boxes of 12 containers	10259



1 oz. Snappies

- 1 oz. container ideal for colostrum collection as milk supply begins to increase
- Can attach directly to a breast pump or be used for hand expression
- Small volume maximizes freezer and storage space

Product	Packaging	Cat No.
1 oz (35 mL)	10 per case	10265
1 oz (35 mL)	100 per case – 10 packs of 10 containers	10256



11 mL Colostrum Collector

- Intended for collection when colostrum first comes in
- 11ml volume minimizes loss of 'precious gold'
- Designed for hand expression

Thermo Scientific™ Snappies™ 11 mL Colostrum Collectors

Product	Packaging	Cat No.
11 mL Colostrum Collector	100 per case	10275

Thermo Scientific™ related products

Ordering information

Product	Description	Packaging	Cat No.
4 oz. (120 mL)	4 oz. lock seal container ideal for small volume mixing applications	250 per case	04HPLS
10 oz. (300 mL)	10 oz. lock seal container best for mixing larger volumes	100 per case	10LSCL
Snappies Snackies	2.5 oz. on-the-go snack containers for mom and child	8 per pack	10261
Snappies Travelers	1.5 oz. containers to conveniently organize small items for home and travel	8 per pack	10260

All Snappies products are proudly made in the USA.

Snappies are BPA free, disposable and recyclable.



Find out more at thermofisher.com/snappies

ThermoFisher
SCIENTIFIC