


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS	
		TUBE	CLOSURE
339650 (BULK)	15	PP	HDPE
339651 (RACKED)	15	PP	HDPE
339652 (BULK)	50	PP	HDPE
339653 (RACKED)	50	PP	HDPE

SIZE.....:	15		50	
NOM. CAP(mL).....:	15		50	
CLOSURE SIZE (mm).....:	17		30	
	INCH	MM	INCH	MM
A (NECK O.D.).....:	0.67	17.0	1.16	29.4
B (NECK I.D.).....:	0.60	15.2	1.09	27.7
C (DIA.).....:	0.62	15.8	1.10	27.9
D (DIA.).....:	0.21	5.3	0.25	6.4
E (HT. @ SHOULDER).....:	3.74	95.0	3.88	98.6
F (TUBE O.A.H.).....:	4.68	118.9	4.50	114.3
G (CLOS DIA.).....:	0.84	21.4	1.34	34.0
H (CLOS. HT.).....:	0.43	11.0	0.51	13.0
J (O.A.H.).....:	4.74	120.4	4.57	116.0

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 33965X	REV. 1111
CK'D	PB	11/30/11			
DESCRIPTION:				CONICAL CENTRIFUGE TUBES	