



Legend:

- 1 Mix Multi mini
06.1250
- 2 Filter housing 10" + prefilter 5 µm/act. carbon
06.5310 + 06.5201
- 3 Filter housing 10" + prefilter 3 µm
06.5310 + 06.5107
- 4 Filter housing 10" + prefilter 0.5 µm
06.5310 + 06.5126
- 5 LabTower EDI 15, 30 or 45
07.4415, 07.4430 or 07.4445

Client / Town:		Toshvin, India		Project:		Scale	
Date		Name		Date		Name	
Drawn: 05.06.07		Opp: Albrecht		Checked: 05.06.07		Date	
Thermo SCIENTIFIC		Thermo Electron LED GmbH		Styckling 3, D-55672 Hiedendorf		Tel. +49 (0)6 021 1 06 99-0	
Water purification system				LabTower EDI			
SZ226				Sheet			
				Sheets			