


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS			
		BOTTLE	CLOSURE	PLUG	O-RING
3140	1000	PC	PPE+HIPS	PP-GLASS FILLED	SILICONE

SIZE CODE.....:	1002/1006/1024	
NOM. CAP.....:	1000	
CLOSURE SIZE (mm).....:	84.5	
	INCH	MM
A (NECK I.D.).....:	2.82 ±.02	71.6 ±0.5
B (HT. W/ CLOSURE).....:	7.69 ±.06	195.3 ±1.5
C (DIA.).....:	3.84 ±.05	97.5 ±1.3
NOM WALL THICKNESS.:	.08	2.0



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 003140-10XX	REV. 0807
CK'D	KOW	8/24/07			

DESCRIPTION: 1L HIGH SPEED CENTRIFUGE BOTTLE