



Rotation, tilting, and vibrating

Versatile and flexible HulaMixer™ Sample Mixer

- Tilt, rotate, and/or vibrate your samples
- Continuous or timed operation
- All settings easily adjustable
- For use at 4°C to 40°C

A lot of labs may not think about it much, but good mixing is important for optimal results. The HulaMixer™ Sample Mixer is perfect for sample preparation with Dynabeads® products and for any other application needing thorough mixing.

Flexible mixing

The multipurpose HulaMixer™ helps mix practically any sample in any tube. Speed is adjustable from 1 to 100 rpm, and it can rotate in three ways: orbital (end-over-end), reciprocal (tilting), and vibrating (vortex). These motions can be performed individually or in



Figure 1. The versatile HulaMixer™ Sample Mixer comes with two tube holder plates to suit all kinds of microtubes and vacutainers (0.5–50 mL).

combinations. The mixer comes with two separate platforms, accommodating tubes ranging from 0.5 to 50 mL.

User-friendly and portable

The HulaMixer™ Sample Mixer is flexible, yet intuitive and easy to operate (Figure 1). You adjust speed, time, and movement combinations to match your changing needs. The HulaMixer™ comes with a manual, a laminated Quick Guide, and AC/DC adaptors. The motor is extremely quiet. The mixer is small; just carry it to the incubator or coldroom.

Personalize your HulaMixer™

Every lab needs a flexible, easy-to-use sample mixer. Or maybe you need a couple of them! Want one all for yourself? The HulaMixer™ comes with a selection of durable stickers, allowing you to add your personal touch to your sample mixer.

Ordering information

Product	Cat. No.
HulaMixer™ Sample Mixer	159-20D

For current prices and further information, visit

www.invitrogen.com/hulamixer.

Stop waiting for an available mixer. Order one for yourself or your group today!