

800/877 Scheduled Customer Maintenance Logsheets

Perform Weekly:	Jan	Feb	Mar	Apr	May	Jun	Jul	Aug	Sep	Oct	Nov	Dec
Replace AC4-1 tube												
Clean origin trenches												
Replace buffer and bleach & clean bottles												
Replace H ₂ O & EtOH/clean external bottles												
Empty waste bottle												
Purge Wells utility												
Perform Monthly*:												
Remove, clean, and oil TC plate												
Remove/inspect TC lid gasket												
Remove/inspect syringes												
Check/fill coolant level												
Perform Every Three Months* and when Replacing TC Plate:												
Replace TC lid gasket												
Perform Every Three Months* (PCR only)												
Replace probe tip												
Replace TC plate and gasket												
Perform Every Six Months* (Sequencing only)												
Replace probe tip												
Replace TC plate and gasket												
*Run comprehensive test after performing maintenance at these intervals												

Instrument S/N:
Other:

Thermal Cycler Sensor Calibration				
37 °C/	95 °C/	Date:	Tool S/N:	Engineer ID:
37 °C/	95 °C/	Date:	Tool S/N:	Engineer ID:
37 °C/	95 °C/	Date:	Tool S/N:	Engineer ID: