



## GeneTitan™ Instrument Scripts for 16 peg HT Array Plates

This ZIP file contains GeneTitan scripts for the 16 peg PM Only HT Array Plates for Human, Mouse and Rat. In order to process the 16 peg PM only HT Array Plates on the GeneTitan Instrument, the files will need to be copied to the following Affymetrix GeneChip Command Console (AGCC) directory

Install directory: **C:\Command\_Console\Fluidics\96F**

The GeneTitan scripts are:

1. GeneTitan™ Instrument script for GeneChip® HT HG-U133+ PM **16 peg** array plate : **550111.protocol**
2. GeneTitan™ Instrument script for GeneChip® HT MG-430 PM **16 peg** array plate : **550112.protocol**
3. GeneTitan™ Instrument script for GeneChip® HT RG-230 PM **16 peg** array plate: **550113.protocol**

After installing the scripts, restart the AGCC GeneTitan Instrument Control software.