

RapID™ STR QUALITY CONTROL REPORT FORM

1. *Record results of biochemical reactions under the expected reactions.
 Corrective Action: See Reverse Page

Date QC performed: _____
 Technologist: _____
 STR lot #: _____
 Expiration Date: _____
 Approved by: _____

Organism	ARG	ESC	MNL	SBL	RAF	INU	GAL	GLU	NAG	PO4	TRY	HPR	LYS	PYR
<i>E. faecalis</i>	+	+	+	+	-	-	-	+	+	(-)	V	-	+	+
ATCC® 29212														
Organism	ARG	ESC	MNL	SBL	RAF	INU	GAL	GLU	NAG	PO4	TRY	HPR	LYS	PYR
<i>E. durans</i>	+	+	-	-	-	-	+	-	+	-	V	+	(-)	+
ATCC® 11576/49479														
Organism	ARG	ESC	MNL	SBL	RAF	INU	GAL	GLU	NAG	PO4	TRY	HPR	LYS	PYR
<i>S. gallolyticus</i> ^a	-	+	+	-	+	+	+	+	V	-	-	-	+	-
ATCC® 9809														
Organism	ARG	ESC	MNL	SBL	RAF	INU	GAL	GLU	NAG	PO4	TRY	HPR	LYS	PYR
<i>S. pyogenes</i> ^a	+	-	-	-	-	-	-	(+)	-	+	+	-	+	+
ATCC® 19615														

+ = Positive V = Variable
 - = Negative (-) = Usually negative
 (+) = Usually positive

*Visual Interpretation Only

^aKey indicator strain

CORRECTIVE ACTION

1. List organism with biochemical reaction that is outside the acceptable criteria:

2. Date corrective action started: _____

3. Corrective action procedure performed: _____

4. Results after corrective action:

Acceptable Reaction

Outside Acceptable Reaction

5. Further corrective action results, if the biochemical reaction is still outside the acceptable criteria:

6. Date corrective results accepted: _____