

**Distribution:** Central File**Date:** 07/10/09**Supersedes:** 18/10/06**OXOID QUALITY ASSURANCE  
PRODUCT SPECIFICATION****KANAMYCIN DISCS****DD0027****Description**

For the presumptive identification of anaerobic gram negative bacteria and particular discrimination of *Bacteroides fragilis*.

**Microbiological Characteristics**

Place discs on to Columbia Blood Agar Base (CM0331) plates enriched with 5% v/v defibrinated horse blood (SR0050), inoculated with the test organisms.

**Reactions after incubation at 37°C for 24-48 hours under anaerobic conditions  
(for details refer to Oxoid Manual - Atmosphere Generation Systems)**

**Zone**

<i>Bacteroides fragilis</i>	ATCC® 25285	Resistant
<i>Fusobacterium varium</i>	ATCC® 8501	Sensitive