

BT-SPEC-0587

Distribution: Central File

Date: 04/08/05

Supersedes: 24/01/05

OXOID QUALITY ASSURANCE

PRODUCT SPECIFICATION

BCET-RPLA (OXOID)

TD950 (20tests)

Description

A kit for the detection of *Bacillus cereus* enterotoxin (diarrhoeal type) by reversed passive latex agglutination.

Components of the kit

TD951	Sensitised latex
TD952	Latex control
TD953	Enterotoxin control (lyophilized)
TD954	Diluent
	Instruction leaflet

For detailed test procedure refer to instruction leaflet X3095.

The product when used as described in the Instruction Leaflet must detect *Bacillus cereus* enterotoxin (diarrhoeal type) at a concentration of 2ng/ml or greater.