

	Document Owner Department: QC	MBD-BT-SPEC-0858
		Page 1 of 1
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LINCOMYCIN MY10 CT0123B		

Product Details

Product Name	LINCOMYCIN
Disc Code	MY10
Disc Concentration	10
CT Code	CT0123B
Shelf Life (years)	3
CLSI Recommended Concentration	Yes
EUCAST Recommended Concentration	No

Microbiological Performance Criteria

CLSI		
Organism	ATCC/NCTC Number	Range (mm)
<i>Staphylococcus aureus</i>	ATCC® 25923	20 - 28
<i>Escherichia coli</i>	ATCC® 25922	0
<i>Pseudomonas aeruginosa</i>	ATCC® 27853	0

References

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests;

- CLSI routine and extended internal quality control for MIC determination and disc diffusion as recommended by CLSI

Other additional challenging strains may be employed

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.

Revision History

Section / Step	Description of Change	Reason for Change	Reference
N/A	Creation of new document	No previous product specification	N/A