	Document Owner Department: QC	BT-SPEC-0762
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
TILMICOSIN TIL15 CT1756B		

Product Details

Product Name	TILMICOSIN
Disc Code	TIL15
Disc Concentration	15µg
CT Code	CT1756B
Shelf Life (years)	1
CLSI Recommended Concentration	Yes
EUCAST Recommended Concentration	No

Microbiological Performance Criteria

CLSI		
Organism	ATCC/NCTC Number	Range (mm)
<i>Staphylococcus aureus</i>	ATCC®25923	17 – 21


References

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests;

- CLSI: VET01 Performance standards for antimicrobial disc and dilution susceptibility tests for bacteria isolated from animals

Other additional challenging strains may be employed.

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.

	Document Owner Department: QC	BT-SPEC-0762
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
TILMICOSIN TIL15 CT1756B		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document/ References	Update to current format. Update of current performance standards	Change control	BT-CC-2802