

	Document Owner Department: QC	MBD-BT-SPEC-0363
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
V FACTOR DISCS DD0004		

V FACTOR DISCS

DD0004

Description

Discs impregnated with Nicotinamide Adenine Dinucleotide for the differentiation of the Haemophilus group of bacteria.

Microbiological Characteristics

Place discs on to Blood Agar Base (CM0055) plates inoculated with the test organisms.

Reactions after incubation at 37°C for 18 hours

<i>Haemophilus influenzae</i>	ATCC®19418	No growth
<i>Haemophilus influenzae</i>	ATCC®49247	No growth
<i>Haemophilus parainfluenzae</i>	ATCC®33392	Growth
<i>Haemophilus haemoglobinophilus</i>	NCTC8540	No growth

	Document Owner Department: QC	MBD-BT-SPEC-0363
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
V FACTOR DISCS DD0004		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document	Update to new document format and correction of minor/typographical errors	N/A	N/A
Microbiological characteristics	Removal of <i>H. influenzae</i> ATCC®9334. Addition of <i>H. haemoglobinophilus</i> NCTC8540.	Change control	MOC-2023-0089