


<b>ThermoFisher</b> SCIENTIFIC	Document Owner Department: QC	MBD-BT-SPEC-0857
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
CARTRIDGE CONTAINER WITH DESICCANT HP0055A DESICCANT SACHETS HP0063A		

## Product Details

<b>Product Name</b>	Thermo Scientific™ Oxoid™ Cartridge Container with Desiccant Thermo Scientific™ Oxoid™ Desiccant Sachets
<b>Product Code</b>	HP0055A HP0063A
<b>Warranty period (years)</b>	N/A
<b>Shelf Life (years)</b>	N/A
<b>Storage</b>	Store single disc cartridges in the cartridge container with desiccant to protect AST discs from moisture. When not in use, store discs in container at 2°C - 8°C. On removal from refrigerator, allow to attain room temperature before use.
<b>Compliance</b>	For use only with Thermo Scientific™ Oxoid™ Antimicrobial Susceptibility Test Discs

	Document Owner Department: QC	MBD-BT-SPEC-0857
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
CARTRIDGE CONTAINER WITH DESICCANT HP0055A DESICCANT SACHETS HP0063A		

## Revision History

Section / Step	Description of Change	Reason for Change	Reference
Product Name	HP0063A added to specification	HP0063A is a component used with HP0055A	IN 1026531
Product Code	HP0063A added to specification	HP0063A is a component used with HP0055A	IN 1026531