




thermo**scientific**



Rapid, reliable results for respiratory season


Results in 15 minutes. Simple procedure. Accurate and reliable. Meets FDA influenza requirements.

ThermoFisher
SCIENTIFIC

Guide proper patient management, increase workflow efficiency

When it comes to Flu A, Flu B and Respiratory Syncytial Virus (RSV) testing, one name stands out: Thermo Scientific™ Xpect™. Results from Xpect Flu A&B and Xpect RSV allow physicians to make a rapid diagnosis with confidence, which is essential for reducing secondary testing, decreasing unnecessary antibiotic use, and guiding appropriate antiviral therapy.

- Easy-to-read results in 15 minutes
- Simple procedure
- Accurate, reliable results
- Long shelf life to adapt to seasonal fluctuations
- Variety of collection and transport options
- Convenient room temperature storage



		Flu A&B	RSV
Performance	Sensitivity	A = 92.2% B = 97.8%	95.6%
	Specificity	100%	94.1%

Annual Analytical Reactivity Testing


On an annual basis, a panel of human influenza viral strains provided by the Centers for Disease Control and Prevention (CDC) is used to evaluate the reactivity of the Xpect Flu A & B Test.

The minimally-reactive concentration is identified as the lowest dilution at which any of the replicates give a positive reading. Testing is performed as per 21 CFR 866.3328 for Influenza virus antigen detection test systems.


Data is published at thermofisher.com/xpectflu.

Simple 2-step, walk-away procedure

Xpect Flu A&B



Xpect RSV



Combine the best in respiratory testing with industry-leading specimen transport

Thermo Scientific™ MicroTest™ M4, M5, M6 and M4RT

MicroTest transport products are proven leaders for the collection and transport of viruses, Chlamydia, Mycoplasma and Ureaplasma. MicroTest products incorporate patented formulations for collection, maintenance, transport and long-term storage.

- The trusted transport media from the media experts backed by over 15 years of proven performance
- Fill volume accommodates multiple tests on the same specimen
- Room temperature storage (select formulations)
- Compatibility with culture, immunoassays, and molecular diagnostics
- Color-coded caps and convenient packaging



Xpect more with one trusted source for all your respiratory testing needs

Ordering information

Product Description	Testing type(s)	Product Number
Xpect		
Xpect Flu A&B	20 Test Kit (Controls included)	R24600
Xpect RSV	30 Test Kit (Controls sold separately)	R24601
Xpect Flu A&B control swab (accessory item for R24600)	20 Flu A+ / B- and 20 Flu A- / B+	R246003
Xpect RSV liquid controls	One each of RSV positive and RSV negative control	R246012
MicroTest		
M4	72 tubes/pack	R12500
M4RT	72 tubes/pack	R12505
M4RT Female/Male Kit	100 kits/pack	R12578
M4 Female/Male Kit	100 kits/pack	R12570
M5	72 tubes/pack	R12515
M6	72 tubes/pack	R12530
MicroTest Micro Tip Flocked Swabs	100 swabs/pack	R12542

Find out more at thermofisher.com/xpectflu