










## Cryo-EM Grant Program

### Accessories List

Here's what you could win when you apply for a Thermo Scientific™ Cryo-EM Grant.

Accessory / Consumable	Part Number	Short Description	Restrictions
 <p>C-Clip Ring (or AutoGrid Ring) (100x)</p>	1036173	Set of 100 rings in which classical 3 mm specimen grids can be mounted. A C-Clip Ring is also referred to as <i>AutoGrid Ring</i> .	Limit 2
 <p>C-Clip (100x)</p>	1036171	Set of 100 clips that secure the specimen grids inside the AutoGrid Rings. An assembly of C-Clip Ring, specimen grid and C-Clip is referred to as <i>AutoGrid</i> .	Limit 4
 <p>C-Clip Insertion Tool</p>	9432 909 97571	Tool to insert a C-Clip into an AutoGrid Ring.	Limit 1

Accessory / Consumable	Part Number	Short Description	Restrictions
<b>AutoGrid Alignment Tool</b> 	9432 909 97641	Tool that assists in aligning the C-Clip Insertion Tool to an AutoGrid Ring.	Limit 1
<b>AutoGrid Container (20x)</b> 	1084591	Set of 20 boxes in which up to 4 AutoGrids can be stored. Each box comes with a lid that can be screwed on to hold the AutoGrids in place and protect them during cryo transfers.	Limit 10
<b>Grid Container Tool</b> 	9432 909 97671	Tool to pick up, transfer, and place both conventional Grid Containers and AutoGrid Containers.	Limit 5
<b>AutoGrid Tweezers</b> 	9432 909 97631	Pair of tweezers that is used to manipulate AutoGrids between an AutoGridbox and a Cassette.	Limit 1
<b>Cassette Tweezers</b> 	9432 909 97651	Pair of tweezers that is used to transfer a Cassette into and out of a Capsule.	Limit 1
<b>Cassette</b> 	1020863	Storage block which can hold up to 12 AutoGrids. The Cassette is the only part that is transferred into and out of the Autoloader.	Limit 1

 Apply now at [thermofisher.com/EM-award](https://thermofisher.com/EM-award)

For research use only. Not for use in diagnostic procedures. For current certifications, visit [thermofisher.com/certifications](https://thermofisher.com/certifications)

© 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. FL0204-EN-12-2022