

Comments for pDEST™17
6354 base pairs

T7 promoter: bases 21-40
Ribosome binding site (RBS): bases 86-92
Initiation ATG: bases 101-103
6xHis tag: bases 113-130
attR1: bases 140-264
Chloramphenicol resistance gene (*Cm^R*): bases 373-1032
ccdB gene: bases 1374-1679
attR2: bases 1720-1844
T7 transcription termination region: bases 1855-1983
bla promoter: bases 2471-2569
Ampicillin resistance (*bla*) gene: bases 2570-3430
pBR322 origin: bases 3575-4248
rop ORF: bases 4619-4810 (c)

(c) = complementary strand

