

Comments for pcDNA™5/FRT
5070 base pairs

CMV promoter: bases 232-819
CMV forward priming site: bases 769-789
T7 promoter/priming site: bases 863-882
Multiple cloning site: bases 895-1010
BGH reverse priming site: bases 1022-1039
BGH polyadenylation sequence: bases 1028-1252
FRT site: bases 1536-1583
Hygromycin resistance gene (no ATG): bases 1591-2611
SV40 early polyadenylation signal: bases 2743-2873
pUC origin: bases 3256-3929 (c)
bla promoter: bases 4935-5033 (c)
Ampicillin resistance (*bla*) gene: bases 4074-4934 (c)

(c) = complementary strand

