

Comments for pcDNA5/FRT/CAT
5858 base pairs

CMV promoter: bases 232-819
 CMV forward priming site: bases 769-789
 T7 promoter/priming site: bases 863-882
 Chloramphenicol acetyl transferase (CAT) gene: bases 1026-1685
 BGH reverse priming site: bases 1810-1827
 BGH polyadenylation signal: bases 1816-2040
 FRT site: bases 2324-2371
 Hygromycin resistance gene (no ATG): bases 2379-3399
 SV40 early polyadenylation signal: bases 3531-3661
 pUC origin: bases 4044-4717 (c)
bla promoter: bases 5723-5821 (c)
 Ampicillin resistance (*bla*) gene: bases 4862-5722 (c)

(c) = complementary strand

