




**Comments for pCMV SPORT-βgal:
7854 nucleotides**

SV40 small T intron and polyA signal: bases 193-555 (complementary strand)

T7 promoter: bases 645-664

lacZ ORF: bases 1009-4149 (complementary strand)

SP6 promoter: bases 4259-4278 (complementary strand)

CMV promoter: bases 4308-4901 (complementary strand)

pUC origin: bases 5390-6063 (complementary strand)

loxP: bases 6115-6148

Ampicillin (*bla*) resistance gene: bases 6250-7110 (complementary strand)

incA: bases 7134-7306

f1 intergenic region: bases 7579-7854 (complementary strand)