




Comments for pDEST22
8930 nucleotides

ADH1 promoter: bases 272-1726
 Nuclear localization signal (NLS): bases 1734-1754
GAL4 DNA activation domain: bases 1761-2105
attR1 site: bases 2121-2245
 Chloramphenicol resistance (Cm^R) gene: bases 2495-3154
ccdB gene: bases 3495-3800
attR2 site: bases 3841-3965
ADH1 transcription termination region: bases 4203-4360
f1 origin: bases 4687-5142
TRP1 gene: bases 5245-5919 (c)
ARS4/CEN6 origin: bases 6455-6972
 Ampicillin (*bla*) resistance gene: bases 7104-7964
 pUC origin: bases 8109-8782