



**Comments for pEXP32/Krev1  
11062 nucleotides**

*ADH1* promoter: bases 103-1557

*GAL4* DNA binding domain: bases 1581-2024

*attB1* site: bases 2037-2061

*Krev1* gene: bases 2082-2630

*attB2* site: bases 2653-2677

*ADH1* transcription termination region: bases 2915-3072

f1 origin: bases 3399-3854

*Leu2* gene: bases 4563-5657

ARS4/CEN6 origin: bases 6385-6903

Gentamicin resistance gene: bases 7248-7781 (c)

pUC origin: bases 8629-9302

Cycloheximide sensitivity (CYH2): bases 10241-10690 (c)

(c) = complementary strand