

Comments for pFastBac™1
4775 base pairs

f1 origin: bases 2-457

Ampicillin resistance (*bla*) gene: bases 589-1449

pUC origin: bases 1594-2267

Tn7R: bases 2511-2735

Gentamicin resistance gene: bases 2802-3335 (c)

Polyhedrin promoter: bases 3904-4032

Multiple cloning site: bases 4037-4142

SV40 polyadenylation signal: bases 4160-4400

Tn7L: bases 4429-4594

(c) = complementary strand

