

pPink-HC Multiple Cloning Site

```

                                5' end of AOX1 mRNA
801 CCCTGTCTTA AACCTTTTTT TTTATCATCA TTATTAGCTT ACTTTCATAA TTGCGACTGG
    5' AOX1 priming site
861 TTCCAATTGA CAAGCTTTTG ATTTTAACGA CTTTAAACGA CAACTTGAGA AGATCAAAAA
                                EcoR I   Rsr II       Sph I       Stu I       Kpn I   Nae I
921 ACAACTAATT ATTCGAAACG GAATTCCGGA CCGGCATGCC AAGGCCTCAG GTACCGGCCG
    Fse I   Swa I   CYC1 transcription termination region
981 GCCATTTAAA TACAGGCCCC TTTTCCTTTG TCGATATCAT GTAATTAGTT ATGTCACGCT

```