

pPink-LC Multiple Cloning Site

```

                                5' end of AOX1 mRNA
                                |
801 CCCTGTCTTA AACCTTTTTT TTTATCATCA TTATTAGCTT ACTTTCATAA TTGCGACTGG
    |
    5' AOX1 priming site
    |
861 TTCCAATTGA CAAGCTTTTG ATTTTAACGA CTTTAAACGA CAACTTGAGA AGATCAAAAA
                                EcoR I   Rsr II       Sph I   Stu I       Kpn I   Nae I   Fse I
                                |       |       |       |       |       |       |
921 ACAACTAATT ATTCGAAACG GAATTCCGGA CCGGCATGCA GGCCTGGTAC CGGCCGGCCA
    |
    Swa I   |
    |       |
981 TTTAAATACA GGCCCTTTT CCTTTGTCGA TATCATGTAA TTAGTTATGT CACGCTTACA
    
```