




**Comments for pYES2.1/V5-His/lacZ:
8964 nucleotides**

GAL1 promoter: bases 1-451

GAL1 Forward priming site: bases 414-437

T7 promoter: bases 475-494

LacZ ORF: bases 522-3578

V5 epitope: bases 3621-3662

V5 C-term Reverse priming site: bases 3630-3650

Polyhistidine (6xHis) region: bases 3672-3689

CYC1 transcription termination signal: bases 3708-3956

pUC origin: bases 4146-4819

Ampicillin resistance gene: bases 4964-5824 (complementary strand)

URA3 promoter: bases 6724-6949 (complementary strand)

URA3 gene: bases 5842-6723 (complementary strand)

2μ origin: bases 6953-8424

f1 origin: bases 8492-8947 (complementary strand)