

**Comments for pYES2.1/V5-His-TOPO™**  
**5886 base pairs**

*GAL1* promoter: bases 1-451

*GAL1* forward priming site: bases 414-437

T7 promoter: bases 475-494

TOPO™ cloning site: bases 512-513

V5 epitope: bases 543-584

V5 C-term reverse priming site: bases 552-572

6xHis tag: bases 594-611

*CYC1* transcription termination signal: bases 633-886

pUC origin: bases 1068-1741

Ampicillin resistance (*bla*) gene: bases 1886-2746 (c)

*URA3* promoter: bases 3647-3872 (c)

*URA3* gene: bases 2764-3645 (c)

2μ origin: bases 3875-5346

f1 origin: bases 5414-5869 (c)

(c) = complementary strand

