

Comments for pYES2
5856 base pairs

GAL1 promoter: bases 1-451

T7 promoter/priming site: bases 475-494

Multiple cloning site: bases 501-600

CYC1 transcription terminator: bases 608-856

pUC origin: bases 1038-1711

Ampicillin resistance (*bla*) gene: bases 1856-2716 (c)

URA3 gene: bases 2734-3841 (c)

2 micron (μ) origin: bases 3845-5316

f1 origin: bases 5384-5839 (c)

(c) = complementary strand

