

CERTIFICATE OF ANALYSIS

Thermo Fisher Scientific's Quality System has been found to conform to Quality Management System Standard ISO9001:2015 by TUV SUD America Inc.
Certificate Number. 951001163

Catalogue Number	J67569
Lot Number	P29L502
Description	Gey's Balanced Salt Solution, with calcium, with magnesium, without phenol
Quality Test/Release Date	03/Apr/2025
Suggested retest date	03/Apr/2027
Country of Origin	Canada

Result Name	Units	Specifications	Test Value
Form		Reported	Clear liquid
pH		Reported	7.4 +/-0.15
Calcium Chloride	mM	Reported	1.5
Disodium Phosphate	mM	Reported	0.8
Glucose	mM	Reported	5.6
KCl	mM	Reported	4.9
Magnesium Chloride	mM	Reported	1.0
Magnesium Sulfate	mM	Reported	0.3
Monopotassium Phosphate	mM	Reported	0.2
Sodium Bicarbonate	mM	Reported	27.0
Sodium Chloride	mM	Reported	121



Seth Christopher
Quality Manager

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third-party data or information associated with the product. Products are for research and development use only. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

Note: The data listed is valid for all package sizes of this lot of this product, expressed as an extension of the catalogue number listed above.