

CERTIFICATE OF ANALYSIS

Thermo Fisher Scientific's Quality System has been found to conform to Quality Management System Standard ISO9001:2015 by Intertek Global
Certificate Number: CERT-0120633

Catalogue Number	Z20
Lot Number	253233
Description	ZINC ACETATE, CERTIFIED
CAS Number	5970-45-6
Quality Test/Release Date	09/Jul/2025
Suggested retest date	03/Jul/2027
Country of Origin	Mexico
BSE/TSE	Above product is free from any substances related to TSE (Transmissible Spongiform Encephalopathy) or BSE (Bovine Spongiform Encephalopathy)

Result Name	Units	Specifications	Test Value
APPEARANCE		REPORT	Clear, white, crystalline material
CALCIUM	%	<= 0.005	0.001
CHLORIDE	%	<= 0.001	<0.001
IDENTIFICATION	PASS/FAIL	= PASS TEST	PASS TEST
INSOLUBLE MATTER	%	<= 0.005	0.001
IRON (Fe)	ppm	<= 5	<5
LEAD (Pb)	%	<= 0.002	<0.002
MAGNESIUM	%	<= 0.005	<0.005
PH 5% SOLUTION @ 25 DEG C		Inclusive Between 6 - 7	6.3
POTASSIUM (K)	%	<= 0.01	<0.01
SODIUM (Na)	%	<= 0.05	<0.05
SULFATE (SO4)	%	<= 0.002	<0.002



Harout Sahagian
QC Manager

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third-party data or information associated with the product. Products are for research use or further manufacturing. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

Note: The data listed is valid for all package sizes of this lot of this product, expressed as an extension of the catalogue number listed above.