

## CERTIFICATE OF ANALYSIS

Thermo Fisher Scientific's Quality System has been found to conform to Quality Management System Standard ISO9001:2015 by TUV SUD America Inc.  
Certificate Number. 951001163

Catalogue Number	J61549
Lot Number	R12L042
Description	Fluorescein disodium salt
CAS Number	518-47-8
Quality Test/Release Date	12/May/2025
Suggested retest date	12/May/2030
Country of Origin	India
Declaration of Origin	Synthetic

Result Name	Units	Specifications	Test Value
Assay (special)	%	≥90.0% (on dry basis)	91.5
Assay (special)	ppm	Residual solvent as isotope alcohol: ≤ 5000 PPM	<320
Assay (special)		Identification A, B, C: Passes tests	Passes Test
Appearance (Color)		Dark brownish-red to brownish-orange, dark orange-red	Brownish
Form		Powder	Powder
Loss on Drying	%	≤ 8.0% (1 hour at 135°C)	0.8
Solution Test		Clear to slightly hazy dark brownish yellow solution with greenish fluorescence (1 mg/ml in water)	Clear dark yellow with greenish fluorescence



John Astarita  
Director, Quality

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third-party data or information associated with the product. Products are for research and development use only. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

Note: The data listed is valid for all package sizes of this lot of this product, expressed as an extension of the catalogue number listed above.