

CERTIFICATE OF ANALYSIS

Thermo Fisher Scientific's Quality System has been found to conform to Quality Management System Standard ISO9001:2015 by Intertek Global Certificate Number. CERT-0120633

Catalogue Number	E138
Lot Number	252206
Description	ETHYL ETHER, ANHYDROUS, ACS GRADE
CAS Number	60-29-7
Quality Test/Release Date	09/May/2025
Suggested retest date	31/Jul/2026
Country of Origin	India
Declaration of Origin	Organic - synthetic
BSE/TSE	This product is manufactured by a synthetic process; no materials of plant or animal origin are incorporated as raw materials, processing aids, or additives.

Result Name	Units	Specifications	Test Value
APPEARANCE		REPORT	Clear, colorless liquid free of suspended matter
ALCOHOL	PASS/FAIL	= PASS TEST	PASS TEST
ASSAY	%	>= 99.0	99.99
CARBONYL COMPOUNDS	%	<= 0.001	<0.001
COLOR	APHA	<= 10	5
EVAPORATION RESIDUE	%	<= 0.001	<0.0001
IDENTIFICATION	PASS/FAIL	= PASS TEST	PASS TEST
PEROXIDE	ppm	<= 1	<1
TITRATABLE ACID	MEQ/G	<= 0.0002	<0.00003
WATER (H2O)	%	<= 0.03	0.002



Matthew Micek
QC Supervisor

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third-party data or information associated with the product. Products are for research use or further manufacturing. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

Note: The data listed is valid for all package sizes of this lot of this product, expressed as an extension of the catalogue number listed above.