

CERTIFICATE OF ANALYSIS

Thermo Fisher Scientific's Quality System has been found to conform to Quality Management System Standard ISO9001:2015 by Intertek. Certificate Number. 2317548

Catalogue Number	41127
Lot Number	A0468371
Description	1-Heptanesulfonic acid sodium salt monohydrate,HPLC grade
CAS Number	207300-90-1
Quality Test/Release Date	06/Jan/2025
Suggested retest date	06/Jan/2030
Country of Origin	INDIA
Declaration of Origin	synthetic

Result Name	Units	Specifications	Test Value
Appearance (Color)		White	White
Appearance (Form)		Powder or flakes	flakes
Infrared spectrum		Conforms	Conforms
Loss on drying	%	7.4 to 9 (100°C, 3 hrs) (vacuum)	7.9 (100°C, 3 hrs) (vacuum)
Titration after Ion exchange	%	>=99.0	99.3
UV		(2.7 g in 100 ml water) pH 3	(2.7 g in 100 ml water) pH 3
UV		at 210 nm A: =<0.1	at 210 nm A:0.0267
UV		at 220 nm A: =<0.05	at 220 nm A:0.0138
UV		at 230 nm A: =<0.03	at 230 nm A:0.0083
UV		at 240 nm A: =<0.02	at 240 nm A:0.005
UV		at 250 nm A: =<0.01	at 250 nm A:0.0043
UV		at 260 nm A: =<0.01	at 260 nm A:0.0058
Solubility		(10 % w/v in water) Clear colorless	(10 % w/v in water) Clear colorless



Geert Torfs
Supervisor, QC

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third-party data or information associated with the product. Products are for research and development use only. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

Note: The data listed is valid for all package sizes of this lot of this product, expressed as an extension of the catalogue number listed above.