



CERTIFICATE OF REGISTRATION

This is to certify that:

Thermo Electron Manufacturing Ltd

Tudor Road Manor Park, Runcorn, Cheshire WA7 1TA UNITED KINGDOM

operates a

QUALITY MANAGEMENT SYSTEM

which complies with the requirements of

ISO 9001:2015

for the following scope

For the design and manufacture of chromatographic media and the supply of chromatographic media and accessories.

Certificate No: QMS42271

Issued: 24 May 2018

Originally Certified: 23 October 2012

Expires: 23 June 2021

Current Certification: 23 May 2018

Nicole Grantham
General Manager SAI Global Certification Services



ISO 9001



WWW.JAS-ANZ.ORG/REGISTER

Registered by:

SAI Global Certification Services Pty Ltd (ACN 108 716 669) 680 George Street Sydney NSW 2000 Australia with SAI Global Assurance Services Ltd Partis House Davy Avenue Knowhill Milton Keynes MK5 8HJ United Kingdom ("SAI Global") and subject to the SAI Global Terms and Conditions for Certification. While all due care and skill was exercised in carrying out this assessment, SAI Global accepts responsibility only for proven negligence. This certificate remains the property of SAI Global and must be returned to SAI Global upon its request. To verify that this certificate is current please refer to SAI Global On-Line Certification register at <http://www.saiglobal.com>

